

# Testfabrics, Inc.

Your source for Quality Consumables



**IMPORTANT**  
*Please Read*  
**Before You Continue**

- Please be sure of the requirements of your test method before ordering.
- If the Multifiber Fabric is to be used for **wash or other agitation testing** (this includes any test using mechanical action, stainless steel balls, and/or a heated soap bath), it is **STRONGLY** recommended that the pieces have a **heat sealed or sewn edge**. This helps to maintain the structural integrity of the test material.
- Cold Cut edges are **ONLY** recommended for stationary testing such as perspiration testing. It is **NOT** advisable to use Cold Cut pieces in any action/agitation testing.
- Under **NO** circumstances should the cutting stripe be removed from the MFF test piece being prepared.
- Please note Testfabrics, Inc. is not responsible for consumer cutting errors and will not accept consumer cut yardage for return. Consumers purchasing full-width MFF yardage assume full responsibility for accurate and appropriate cutting, heat sealing, or sewing of test pieces.
- Pre-cut pieces save both time and labor helping to keep the priority on the test and not the materials required. Pre-cut heat-sealed, sewn edge, and cold cut pieces are all available in a variety of pack sizes for convenient ordering.

---

Testfabrics, Inc. | 415 Delaware Ave, West Pittston, PA 18643 USA  
Tel: +1 (570) 603 0432 | Fax: +1 (570) 603 0433 | Email: [info@testfabrics.com](mailto:info@testfabrics.com)

[www.testfabrics.com](http://www.testfabrics.com)